**SERAFIT®**

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POLYGLYCOLIC ACID SUTURE MATERIAL**• SERAFIT®****DESCRIPTION**

SERAFIT® is a braided, synthetic, absorbable, sterile surgical suture material. It is made from the homopolymer of glycolic acid and has the molecular formula $[C_2O_2H_2]_n$. It is supplied dyed with D+C violet no. 2, Colour Index No. 60725 or undyed. The suture material satisfies the essential requirements of Directive 93/42/EEC on medical devices and the harmonized standard of the European Pharmacopoeia for "sterile, braided, absorbable, synthetic threads (fila resorbilia synthetica torta sterilia)".

• SERAFIT®

polyglycolic acid, braided, coated with calcium stearate and polycaprolactone or polyglycolic acid/poly lactide

INDICATIONS

SERAFIT® is intended for use in adaptation of soft tissues and in ligatures where medium-term absorbable thread material is required, including ophthalmic surgery and microsurgery. Inadequate data are available as yet on use in the central vascular system.

USE

In the choice and use of suture material, due account should be taken of the condition of the patient, the experience of the doctor, the surgical technique, and the type and size of the wound. Standard knot-tying techniques should be used to ensure knot-holding security. In order to prevent needle damage during suturing it is recommended that the needle be grasped at an area between the end of its first third (as seen from its end) and its midpoint.

ACTION

The function of a medical suture in wound care is

to attach tissue to tissue. Ligatures are a special type of suture intended to close hollow organs. SERAFIT® is degraded by hydrolysis in tissues to form glycolic acid, which is then metabolized in the body. This absorption becomes manifest firstly as a reduction in tensile strength and later as a loss of mass of the thread. The breaking load falls to 60–70% of its initial value after 14 days and to 25–45% of its initial value after 21 days. Degradation is complete after 60 to 90 days.

CONTRAINDICATIONS

SERAFIT® must not be used where more prolonged or long-term stability of the suture material is required.

WARNINGS / PRECAUTIONS / INTERACTIONS

The requirements of the suture material and the techniques required vary with the indication. Therefore, the user should be familiar with the surgical techniques to be used. In the choice of suture material, account should be taken of the in-vivo characteristics of the suture material concerned. In patients with delayed wound healing SERAFIT® should be used only after careful consideration. Like all foreign bodies, SERAFIT® can give rise to calculus formation when in prolonged contact with salt-containing solutions (e.g. in the biliary or urinary tract). Additional use of nonabsorbable suture materials or suitable immobilization techniques may be required with sutures that are under tension or require greater support. Use of suitable ancillary measures is required with sutures in infected tissue, as the multifilament structure of the thread can prolong infection. Because of the increased risk of infection associated with the multifilament structure of the thread, skin sutures should if possible be intracutaneous and the suture material, in particular the knots, should be placed as deep as possible.

Surgical instruments should be handled with care in order to avoid damage to the thread.

Bending of needles or grasping them outside of the recommended area can result in needle breakage. Because of the risk of infection from needle

pricks, needles should be handled and disposed of with particular care.

ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of SERAFIT®: transient local irritation, transient inflammatory reactions to the foreign body; rarely, formation of suture fistulas or granulomas, formation of calculus after prolonged contact with salt-containing solutions.

STERILITY

SERAFIT® is supplied sterile (sterilization method: ethylene oxide). Do not resterilize, critical changes in suitability for use possible (more details available on request)! Do not use damaged packages. Packages that have been opened for use but then not used are to be discarded.


STORAGE CONDITIONS

SERAFIT® should be stored at not more than 25°C, in a clean and dry place. Do not use after the expiry date.

PRESENTATION

Suture materials of SERAFIT® are supplied in various gauge sizes and lengths, dyed and undyed and with various types of needle or without needles. The threads or combinations are supplied in individual packs and multiple packs in a sterile package. More precise details are given in the catalogue.

EXPLANATION OF SYMBOLS ON THE PACKAGE

 **CE**₁₀₁₄ CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

Surgical suture, absorbable, braided, coated, undyed

Surgical suture, absorbable, braided, coated, dyed

Surgical suture, absorbable, braided, stiffened, dyed

Polyglycolic acid

Antireflective needle

Detachable needle

Needle identification symbols and abbreviations are explained in the catalogue.

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SERAG-WIESSNER GmbH & Co. KG