SERAPREN®

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NAHTMATERIAL aus
POLYPROPYLEN

LT VARTOJIMO INSTRUKCIJA

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POLYPROPYLENE SUTURE MATERIAL • SERAPREN® blue

DESCRIPTION

SERAPREN® is a synthetic, monofilament, nonabsorbable, sterile surgical suture. The raw material for the suture is composed of isotactic polypropylene with the empirical formula $[C_3H_6]_n$. SERAPREN® blue is dyed with phthalocyanine blue, Colour Index Nr. 74160 or Solvent Blue, Colour Index No. 61568.

Propylene sutures are available in combination with swaged high-grade stainless steel needles.

They are also available in different variations with special accessories for intracuticular suturing. The suture material satisfies the essential requirements of the Medical Device Directive 93/42/EEC and apart from marginally deviations with respect to the diameter (suture material < EP 1) the harmonised standard of the European Pharmacopoeia "Sterile nonabsorbable sutures (fila non resorbilia sterilia)", and of the monograph "Nonabsorbable surgical suture" of the latest edition of the United States Pharmacopoeia (USP).

INDICATIONS

SERAPREN® is indicated for use in soft-tissue approximation and/or for ligation in cases where nonabsorbable suture material is indicated, including use in cardiovascular and neurological procedures as well as in microsurgery and ophthalmic procedures. Polypropylene sutures can also be used as stay sutures and for tagging.

USE

Due account should be taken of the condition of the patient, the experience of the surgeon, the surgical technique, and the size of the wound when selecting and using the suture material. Standard knot-tying techniques should be used to ensure knot security. In order to prevent damage to the needle while suturing it is recommended that the needle be grasped towards the swaged end at a point between the end of the first third and the mid-point.

ACTION

The function of the surgical suture in wound care is to create tissue-to-tissue approximation or to attach tissue to allogeneic material. The ligature, or tie, as a special form of suture serves to occlude lumens and hollow organs.

After an initial minimal tissue reaction, polypropylene suture material is slowly and progressively encapsulated by connective tissue. It is nonabsorbable and therefore retains its tensile strength indefinitely. Due to its monofilament structure SERAPREN® may also be used in tissue at risk of infection.

CONTRAINDICATIONS

None known.

WARNINGS / PRECAUTIONS / INTERACTIONS

The demands placed on the suture material and the required techniques will vary with the indication. The user should therefore be familiar with the various surgical techniques and procedures and should consider the in-vivo characteristics when selecting the suture material. As with any foreign body SERAPREN® may give rise to calculus formation with prolonged contact with the biliary or urinary tracts. The use of suitable ancillary measures is recommended when suturing within infected tissue. The surgical instruments should be handled with care to avoid damage to the suture. Bending the needles or grasping them outside the recommended area may result in needle breakage. Needles should be handled and disposed of with particular care due to the risk of infection from needle stick injuries.

ADVERSE REACTIONS

Even with correct use of SERAPREN® the following adverse effects may occur: Transient localized irritation, transient inflammatory foreign-body reactions, the formation of fistula or granuloma.

STERILITY

SERAPREN® is supplied sterile (sterilization method: ethylene oxide). Do not resterilize, critical changes in suitability for use possible (more details available on request)! Do not use damaged packs! Discard unused opened packs!

STORAGE CONDITIONS

SERAPREN® should be stored in a clean dry place at room temperature. Do not use after expiry date.

INSTRUCTIONS FOR REMOVAL FROM THE PACK

Minibox: For trouble-free removal the suture should be withdrawn from the minibox vertically and not at an angle.

Minibox – microsutures in protective tube: Detach the tube from the minibox fixation; grasp

needle with needle-holder and withdraw the suture from the protective tube; if a **shorter suture** is needed, it should be cut to the required length together with the protective tube; with **double-armed sutures** 2 separate sutures can be removed, if required, by cutting the suture loop at the end of the protective tube.


PRESENTATION

SERAPREN® is available in a variety of strengths and lengths in the form of needle-suture combinations attached to various needle types. Pledgets to reinforce the body's own structures may be enclosed with SERAPREN®. The combinations are supplied sterile in single and multiple packs (multipack). The intracuticular suture is provided together with special accessories. Please see catalogue for precise details.

EXPLANATION OF SYMBOLS ON THE PACKAGE

 1275 CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

 Surgical suture, nonabsorbable, monofilament, dyed

 Polypropylene

 Antireflective needle

 Detachable needle

 Pledget (length times width times thickness [mm])

Needle identification symbols and abbreviations are explained in the catalogue.

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SERAG-WIESSNER GmbH & Co. KG