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NAHTMATERIAL aus POLYPROPYLEN

**en INSTRUCTIONS FOR USE**

POLYPROPYLENE SUTURE MATERIAL

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MATERIEL DE SUTURE en POLYPROPYLÈNE

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MATERIALE DA SUTURA in POLIPROPILENE

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MATERIAL DE SUTURA de POLIPROPILENO

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MATERIAL DE SUTURA à base de POLIPROPILENO

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## ● SERAPREN®

**WARNINGS AND PRECAUTIONS**

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified personnel, that is responsible to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: ethylene oxide).

- The product is intended for single use only.

- Do not resterilize or re-use the product as critical changes in performance and safety are possible.

- Do not use after expiration date.

- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or unintentionally opened prior to use.

- Use the product only for the intended indication.

**DESCRIPTION AND INTENDED PURPOSE**

SERAPREN® is a synthetic, monofilament, non-absorbable surgical suture. The raw material for the suture is composed of isotactic polypropylene. SERAPREN® is dyed with phthalocyanine blue, Colour Index No. 74160. It is available in various gauge sizes and lengths.

Polypropylene sutures are supplied in combination with various types of needles made of surgical steel. Pledgets to reinforce the body's own structures may be enclosed with SERAPREN®.

The suture material complies with the current edition of the European Pharmacopoeia monograph (EP) "Sutures, sterile non-absorbable (fila non resorbilia sterilia)".

SERAPREN® is indicated for use in soft-tissue approximation and/or for ligation in cases where non-absorbable suture material is indicated, including use in cardiovascular and neurological procedures as well as in microsurgery and ophthalmic procedures. Polypropylene sutures can also be used as stay sutures and for tagging.

**INSTRUCTIONS FOR USE**

In the choice and use of suture material, account should be taken of the condition of the patient, the experience of the qualified personnel, the surgical technique, and the type

and size of the wound. Standard knot-tying techniques should be used to ensure knot security. In order to prevent damage to the needle during suturing it is recommended that the needle be grasped at an area between the end of its first third (as seen from its end) and its midpoint.

### **ACTION**

The function of a medical suture in wound care is to attach tissue to tissue or tissue to foreign material. Ligatures are a special type of suture intended to close hollow organs. After an initial minimal tissue reaction, polypropylene suture material is slowly and progressively encapsulated by connective tissue. Due to its monofilament structure SERAPREN® may also be used in tissue at risk of infection.

### **CONTRAINDICATIONS**

The device must not be used in patients with allergies to any of the product components.

### **RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS**

The requirements of the suture material and the techniques required vary with the indication. In the choice of suture material, account should be taken of the characteristics of the suture material concerned. Like all foreign bodies, the device can give rise to calculus formation when in pro-longed contact with salt-containing solutions (e.g. in the biliary or urinary tract). The use of suitable ancillary measures is recommended when suturing within infected tissue. The surgical instruments should be handled with care to avoid damage to the suture. Bending the needles or grasping them outside the recommended area may result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled with particular care.

### **ADVERSE REACTIONS**

Even with correct use of SERAPREN® the following adverse effects may occur: localized irritation, inflammatory reactions to the foreign body, formation of suture fistulas or granulomas.

Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

### **STORAGE CONDITIONS**

The storage conditions apply as indicated on the label.

### **DISPOSAL INFORMATION**

Disposal of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

### **INSTRUCTIONS FOR REMOVAL FROM THE PACK**

**Minibox:** For trouble-free removal the suture should be withdrawn from the minibox vertically and not at an angle.

**Minibox – microsutures in protective tube:** Detach the tube from the minibox fixation; grasp needle with needle-holder and withdraw the suture from the protective tube; if a **shorter suture** is needed, it should be cut to the required length together with the protective tube; with **double-armed sutures** 2 separate sutures can be removed, if required, by cutting the suture loop at the end of the protective tube.

### **MEANING OF THE SYMBOL**

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG