**SERANOX®****de** GEBRAUCHSANWEISUNGNAHTMATERIAL aus
STAHL DRAHT**pl** INSTRUKCJA OBSŁUGINICI CHIRURGICZNE z
DRUTU STALOWEGO**en** INSTRUCTIONS FOR USESUTURE MATERIAL of
STEEL WIRE**ru** ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮШОВНЫЙ МАТЕРИАЛ из
СТАЛЬНОЙ ПРОВОЛОКИ**fr** MODE D'EMPLOIMATERIAU DE SUTURE en
FIL D'ACIER**nl** GEBRUIKSAANWIJZINGCHIRURGISCH HECHTMATERIAAL uit
STAALDRAAD**it** ISTRUZIONI PER L'USOMATERIALE DA SUTURA in
FILO D'ACCIAIO**ro** INSTRUCȚIUNI DE UTILIZAREMATERIALE DE SUTURĂ din
SĂRMĂ DE OȚEL**es** INSTRUCCIONES DE USOMATERIAL DE SUTURA de
ALAMBRE DE ACERO**lv** LIETOŠANAS INSTRUKCIJAŠŪJAMAIS MATERIĀLS no
TĒRAUDA STIEPLES**cs** NÁVOD K POUŽITÍCHIRURGICKÝ ŠICÍ MATERIÁL z
OCELOVÉHO DRÁTU**mk** УПАТСТВО ЗА УПОТРЕБАМАТЕРИЈАЛИ ЗА ШИЕЊЕ од
ЧЕЛИЧНА ЖИЦА**tr** KULLANIM TALİMATIÇELİK TELDEN
DİKİŞ MATERYALLERİ**lt** VARTOJIMO INSTRUKCIJASIŪLAI iš
PLIENO VIELOS

NAHTMATERIAL aus STAHLDRAHT

- SERANOX® gezwirnt
- SERANOX® gezwirnt, ummantelt
- SERANOX® monofil
- SERANOX® Ti monofil

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DESCRIPTION

SERANOX® is a monofilament or twisted nonabsorbable sterile surgical suture material. The basic material is stainless steel of quality 316L. The twisted version of SERANOX® can be coated with polyethylene, while the monofilament variant can be coated with titanium nitride.

Our steel wire suture materials are supplied as combinations with stainless steel needles or without needles. They are also supplied in a number of variants with specific accessories for retention sutures and tenorrhaphy. They satisfy the essential requirements of Directive 93/42/EEC on medical devices and the harmonized standard of the European Pharmacopoeia (EP) for "suture, sterile non-absorbable (fila non resorbilia sterilia)" with the exception of the chemical composition that corresponds to the United States Pharmacopoeia monograph "Nonabsorbable surgical suture". Twisted, coated SERANOX® complies with the specifications of the EP for the diameter without a coating; the size with a coating can exceed this and is given in mm.

INDICATIONS

Suture material of steel wire is suitable for adaptation of tissues in which nonabsorbable suture material is indicated and high tensile strength is required, in particular bone, tendons, and wounds that are under great tension.

USE

In the choice and use of suture material, due account should be taken of the condition of the patient, the experience of the doctor, the surgical technique, and the size of the wound. Either standard knot-tying techniques or twisting the

ends of the suture can be used to ensure knot-holding security. In order to prevent needle damage during suturing it is recommended that the needle be grasped at a point between the end of its first third (as seen from its end) and its midpoint. When applying a tendon suture according to Lengemann, a separation of 4-5 mm must be maintained between the fixation clip and the skin or the supporting disk, as suture material pulls into the clip housing upon fixation.

ACTION

The function of a medical suture in wound care is to attach tissue to tissue or tissue to allogeneic material.

After a slight initial tissue reaction, suture material of steel wire is slowly and progressively encapsulated by connective tissue. It is nonabsorbable and retains its tensile strength indefinitely. Because of its thread structure, monofilament SERANOX® can also be used in tissue at risk of infection.

CONTRAINDICATIONS

Do not use stainless steel suture material for patients who are allergic to stainless steel of quality 316L (nickel, chrome, etc.).

WARNINGS / PRECAUTIONS / INTERACTIONS

The requirements of the suture material and the techniques needed vary with the indication. Therefore, the user should be familiar with the surgical techniques to be employed. In the choice of suture material, account should be taken of the in-vivo characteristics of the suture material concerned. Due to the higher stiffness the material may be more difficult to knot; throwing large loops will enhance the knotting and reduce the probability of breaking. Care should be taken to ensure that the ends of the sutures are placed properly to avoid irritation of the tissue. When monofilament steel wire sutures are used in infected tissue the use of suitable ancillary measures is recommended; multifilament SERANOX® should not be used in infected tissue, as its thread structure can prolong infection. Surgical instru-

ments should be handled with care in order to avoid damage to threads. Bending of needles or grasping them outside of the recommended area can result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled and disposed of with particular care. Contact corrosion may occur with SERANOX® in combination with other metals. Furthermore, the steel may cause anomalies in certain imaging techniques such as MRI.

ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of SERANOX®: transient local irritation, inflammatory reactions to the foreign body; formation of suture fistulas or granulomas.

STERILITY

SERANOX® is supplied sterile (sterilization procedure: ethylene oxide gas). Do not resterilize, critical changes in suitability for use possible (more details available on request)! Discard unused opened packs!

STORAGE CONDITIONS

The suture material should be kept dry and clean at room temperature. Do not use after the expiry date!

PRESENTATION

SERANOX® is supplied in various thread structures, diameters, and lengths and with various types of needle or without needles. The threads and combinations are supplied in individual packs and multiple packs in a sterile package, and in the case of retention sutures and tenorrhaphy can be supplied with specific accessories. More precise details are given in the catalogue.

EXPLANATION OF SYMBOLS ON THE PACKAGE

CE 1014 CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

Surgical suture, nonabsorbable, monofilament, undyed

Surgical suture, nonabsorbable, monofilament, coated, undyed

Surgical suture, nonabsorbable, twisted, undyed

Surgical suture, nonabsorbable, twisted, coated, undyed

Steelwire

Needle identification symbols and abbreviations are explained in the catalogue.

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SERAG-WIESSNER GmbH & Co. KG