SERANOX®

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NAHTMATERIAL aus STAHLDRAHT

en INSTRUCTIONS FOR USE

SUTURE MATERIAL of STEEL WIRE

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MATERIAU DE SUTURE en FIL D'ACIER

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- SERANOX[®]
- SERANOX[®] Ti

WARNINGS AND PRECAUTIONS

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified personnel, that is responsible to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: ethylene oxide).

- The product is intended for single use only.

- Do not resterilize or re-use the product as critical changes in performance and safety are possible.

- Do not use after expiration date.

- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or unintentionally opened prior to use.

- Use the product only for the intended indication.

DESCRIPTION AND INTENDED PURPOSE

SERANOX[®] is a monofilament non-absorbable surgical suture material. The basic material is stainless steel of quality 1.4404. SERANOX[®] can be coated with titanium nitride.

Steel wire suture materials are available in various gauge sizes and length. They are supplied with various types of needles made of surgical steel or without needles. The suture material fulfills the physical properties of the current edition of the European Pharmacopoeia monograph (EP) for "Sutures, sterile non-absorbable (Fila non resorbilia sterilia)".

Suture material of steel wire is suitable for adaptation of tissues in which non-absorbable suture material is indicated and high tensile strength is required.

INSTRUCTIONS FOR USE

In the choice and use of suture material, due account should be taken of the condition of the patient, the experience of the qualified personnel, the surgical procedure, and the type and size of the wound. Either standard knot-tying techniques or twisting the ends of the suture can be used to ensure knot-holding security. In order to prevent

needle damage during suturing it is recommended that the needle be grasped at a point between the end of its first third (as seen from its end) and its midpoint.

ACTION

After a slight initial tissue reaction, suture material of steel wire is slowly and progressively encapsulated by connective tissue. It is non-absorbable and retains its tensile strength indefinitely.

CONTRAINDICATIONS

The device is not intended specifically for use in direct contact with the central nervous system, the heart or the central circulatory system. The device must not be used in patients with allergies to any of the product components.

RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS

The requirements of the suture material and the techniques needed vary with the indication. In the choice of suture material, account should be taken of the characteristics of the suture material concerned. Like all foreign bodies, the device can give rise to calculus formation when in prolonged contact with salt-containing solutions (e.g. in the biliary or urinary tract). Due to the higher stiffness the material may be more difficult to knot; throwing large loops will enhance the knotting and reduce the probability of breaking. Care should be taken to ensure that the ends of the sutures are placed properly to avoid irritation of the tissue. Because of its thread structure, SERANOX[®] can also be used in tissue at risk of infection. In the case of sutures in infected tissue the use of suitable ancillary measures is called for. Surgical instruments should be handled with care in order to avoid damage to threads. Bending of needles or grasping them outside of the recommended area can result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled with particular care.

Contact corrosion may occur with SERANOX[®] in combination with other metals.

Furthermore, the steel may cause anomalies in certain imaging techniques such as MRI. Depending on the type of device used and the field intensity and duration, local heating, image artifacts and dislocation may occur during an MRI procedure.

ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of the device: local irritation, inflammatory reactions, formation of foreign-body giant cells, suture fistulas, or granulomas.

Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

STORAGE CONDITIONS

The storage conditions apply as indicated on the label.

DISPOSAL INFORMATION

Disposal of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

MEANING OF THE SYMBOLS

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG