**de GEBRAUCHSANWEISUNG****Synthetisches, steriles, chirurgisches, textiles Implantat****en INSTRUCTIONS FOR USE****Synthetic sterile surgical textile implant****fr MODE D'EMPLOI****Implant chirurgical textile synthétique stérile****it ISTRUZIONI PER L'USO****Impianto tessile chirurgico sintetico sterile****es INSTRUCCIONES DE USO****Implante quirúrgico textil sintético estéril****pt INSTRUÇÕES DE UTILIZAÇÃO****Implante têxtil cirúrgico estéril sintético****pl INSTRUKCJA UŻYCIA****Syntetyczny jałowy tekstylny implant chirurgiczny****cs NÁVOD K POUŽITÍ****Syntetický sterilní chirurgický textilní implantát****It NAUDOJIMO INSTRUKCIJA****Sintetinis sterilus chirurginis tekstilinis implantas****Iv LIETOŠANAS INSTRUKCIJA****Sintētisks, sterils, ķirurģisks auduma implants****ro INSTRUCȚIUNI DE UTILIZARE****Implant chirurgical steril din material textil sintetic****mk УПАТСТВО ЗА УПОТРЕБА****Синтетички стерилен хируршки текстилен имплант**

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## Synthetic, sterile, surgical, textile implant

### WARNINGS AND PRECAUTIONS

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified and trained personnel. It is the physician's responsibility to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: ethylene oxide). 
- The product is intended for single use only. 
- Do not re-sterilize or reuse the product as critical changes in performance and safety are possible. 
- Do not use after expiration date. 
- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or accidentally opened prior to use. 
- Use the product only for the intended indication.
- Please contact the manufacturer for further information.

### DESCRIPTION AND INTENDED PURPOSE

SERAG-WIESSNER textile implants are surgical meshes and tapes made of synthetically produced polypropylene threads. These basic threads can be multifilament, monofilament or as partially absorbable bicomponent threads with absorbable coating made of glycolide and  $\epsilon$ -caprolactone. Surgical suture materials can be integrated as supporting threads. These are not intended to remain in the body. The textile implants can be dyed with D+C violet No. 2, Colour Index No. 60725, phthalocyanine blue, Colour Index No. 74160 or Pigment Black 11, Colour Index No. 77499 (iron oxide) over the complete surface or on stripes for orientation. The dyeing with iron oxide enables a representation in the MRI. The implant does not affect magnetic resonance (MR) environments and is not affected in its functionality.

Textile implants are indicated for the surgical support, bridging and/or reinforcement of the body's own structures.

The areas of application concern hernia surgery (SERAMESH®), urology and urogynecology (SERASIS®, SERATOM®, SERATIM® and SERATEX®) as well as breast reconstruction in senology (SERAGYN®).

Different shapes and structures are available for different applications. The respective basic structure can be identified by the letter combinations of the product designation. On request, we can provide you with special data sheets for the individual versions.

Both product labels with documentation-relevant data for the patient surgery report and an implant passport for the patient are included with every textile implant.

### INSTRUCTIONS FOR USE

Users should be familiar with surgical techniques and risk management for complications. Studies have shown that the main risk for erosion or other complications (e.g. bleeding, organ and vascular injuries or recurrences) can lie in the intraoperative

procedure, which is why we offer workshops. Please contact us.

Particular attention must be paid to the correct positioning of the textile implant. If the position is incorrect, correction or removal and correct insertion with a new textile implant is necessary, as otherwise late and further damage (e.g. erosion of the mucous membrane) is possible.

If the textile implant is to be cut individually, it must be tested whether the integrity of the mesh has been impaired.

The implantation should be performed with such care that damage to larger blood vessels, nerves and organs is avoided. The consideration of local anatomy, individual patient status and correct instrument passage minimizes risks.

Special surgical instruments from SERAG-WIESSNER are available for the implantation of textile implants. Please refer to the separate product information.

Surgical instruments should be handled carefully to avoid damaging the textile implant and any sutures used for fixation. The tapes, arms or supporting threads on the textile implants must be inserted carefully into the instrument eyelet so that there is no threading of single thread systems, which can lead to a torn of the band or arm. Textile implants can be mechanically damaged, for example, by clamps, clips or clamps. The textile implant must not be attached with clips or sutures in the outermost mesh at the implant margin. A distance of at least 6 mm or 2 pores from the edge of the mesh must be maintained.

#### Urology and Urogynaecology

It must be ensured that the textile implant has been placed with the correct tension (usually tension-free, but with SERATIM® under tension) and without torsion/twisting. After an overcorrection, a temporary or permanent stenosis of the urinary tract as well as functional disturbances of the bowel and bladder voiding can occur.

Adjustment of the tape tension is possible with the SERASIS® and SERATIM® incontinence implants via the tape ends or special supporting threads. The supporting threads are not intended to remain in the body.

For surgery around the bladder, cystoscopy is recommended to confirm the integrity of the bladder or to detect a bladder perforation or misplacement of the implant. This applies similarly to the intestine or other hollow organs.

#### Reconstructive breast reconstruction

The SERAGYN® insert should be crease-free.

### **PATIENT INFORMATION**

Before the surgery, the patient must be informed individually by a physician about the risks involved in the implantation of surgical meshes and/or tapes. This also includes the contraindications and undesirable side effects mentioned in these instructions for use. Furthermore, adequate postoperative behaviour must be explained by the physician. The general, nationally valid recommendations, e.g. of the medical professional associations, must be observed.

In the event of postoperative complications, the patient must be informed to consult a physician. It should be noted that even years after the operation mesh erosions or mesh migrations or adhesions may occur.

#### Hernia surgery

After hernia surgery, the patient must not lift heavy objects or engage in strenuous sports (a particular swelling at or near the wound site or implant may be a sign of a recurrent hernia).

## Urology and Urogynaecology

After the operation in the pelvis and depending on the complexity of the operation, the patient must not perform any heavy lifting, strenuous sports or sexual intercourse.

## Breast reconstruction

After breast reconstruction with SERAGYN® it is recommended to wear a bra which supports the breast (e.g. Stuttgart-Belt). Heavy lifting, physical work and sporting activities should be avoided.

## CONTRAINDICATIONS

### General

The general contraindications for surgical treatment with textile implants apply.

When using these implants on young patients whose body growth has not yet been completed, or on patients with present or planned pregnancies, it should be taken into account that the textile implants have only limited elasticity and do not grow with the patient.

In patients with an impaired immune system, multiple allergies or hypersensitivity, especially if allergy to the implant material is known, the implantation must be carefully considered. An allergy test should be performed beforehand.

In patients with haemorrhagic diathesis (increased tendency to bleed), the textile implant should only be used after a detailed risk-benefit consideration or appropriate measures have been taken.

Textile implants should be used with caution in patients with delayed wound healing (e.g. diabetics).

Use in infected tissues may promote the spread of the infection.

## Hernia surgery

The lightweight SERAMESH® XO should not be used to repair structures where high pressures can act on the mesh, such as larger abdominal wall hernias.

## Urology and Urogynaecology

Do not use the urological or gynecological pelvic floor implant in patients with urinary tract infections or relative obstruction of the upper urinary tract.

## Reconstructive breast reconstruction

Radiotherapy after mesh insertion is contraindicated as it may lead to reduced neovascularization and mesh integration.

## RESTRICTIONS AND LIMITATIONS OF USE

This product is not intended for use in direct contact with the heart, central circulatory system or the central nervous system.

## UNDESIRABLE SIDE EFFECTS

### General

The following adverse effects may occur:

- Thrombosis/embolism during or after the operation
- Critical bleeding during or after surgery
- Damage (punctures or lacerations) to blood vessels, nerves, structures or organs, including bladder, ureter or intestines.
- Haematoma and seroma formation

- Swelling and/or with temporary local irritation in the surrounding area of the wound site / implant
- Foreign body reaction (e.g. granuloma formation), hypersensitivity, allergic reaction to the implant or chronic inflammation of the tissue. These can cause fibrosis.
- Rejection reactions to the textile implant
- Infections (e.g. wound infection, abscess, tissue necrosis, sepsis)
- Excessive scarring. They may need a new surgical intervention.
- Mesh erosion or migration (e.g. perforation or erosion of tissue, structures, organs with sometimes serious impacts)
- Persistent discomfort after implantation of the textile implant (e.g. chronic pain). These can lead to psychological stress (e.g. depression).
- Serious complications leading to the death of the patient

In addition, depending on the area of application of the textile implants, the following adverse effects must be taken into account:

#### Hernia surgery

- Fistula formation and adhesions
- Acute and/or chronic pain (e.g. neuromuscular problems in the groin, thigh, leg, abdomen), numbness, post-herniorrhaphic neuralgia
- Intestinal obstruction
- Recurrent hernia

#### Urology and Urogynaecology

- Fistula formation
- Acute and/or chronic pain (e.g. neuromuscular problems in the groin, thigh, leg, abdomen), numbness
- Dyspareunia (pain during sexual intercourse)
- (de novo) functional disorders of bowel and bladder voiding (e.g. overactive bladder, urinary retention, constipation), recurrent incontinence
- Atypical vaginal discharge
- Recurrent prolapse (e.g. vaginal prolapse, rectal prolapse, uterine prolapse)
- Osteomyelitis

#### Reconstructive breast reconstruction

- Capsule formation and/or contracture (e.g. foreign body reaction)
- Acute and/or chronic pain (e.g. nipples or neuromuscular problems in the chest area), sensitivity impairment in the chest area/on the nipples

Before releasing the patient from the hospital, check for symptoms and signs of infection, immune reactions (e.g. allergic reaction), bleeding and impairment of general condition.

The treatment of complications may require another surgical procedure, during which the implant may have to be (partially) removed (e.g. severe infections).

Users and patients are encouraged to report any serious incidents that may have occurred in connection with a SERAG-WIESSNER textile implant to the manufacturer and/or the competent authority.

**STORAGE CONDITIONS**

The storage conditions apply as indicated on the label.

**DISPOSAL**

Disposal of used/unused products or waste material should be handled with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

**PRESENTATION**

The textile implants from SERAG-WIESSNER are supplied in various designs and sizes. If required, they are available with incisions for adaptation to anatomical structures. The dimensions < 15 x 15 cm are supplied flat, the rest are packed folded. Please refer to our catalogue for more detailed information.

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SERAG-WIESSNER GmbH & Co. KG