# KNOCHENWACHS

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## BONE WAX

## DESCRIPTION

Bone wax is a sterile, non-resorbable, soft pliable mass with the following composition:

white beeswax BP	89 % w/w
isopropyl palmitate BP	11 % w/w

Bone wax is supplied sterile in rectangular slabs.

## INDICATIONS

Bone wax is designed to control local bleeding from raw bone surfaces.

## USE

Bone wax is applied sparingly to the affected bony structures. Kneading the wax under sterile conditions will warm it until the desired consistency is obtained.

## ACTION

Bone wax prevents bleeding from raw bone surfaces by forming a mechanical barrier.

## CONTRAINDICATIONS

Bone wax should not be used in the presence of likelihood of infection in the area to be treated or where rapid bone regeneration is desired.

## WARNINGS / PRECAUTIONS / INTERACTIONS

Bone wax should be prepared and used under strictly aseptic conditions. It should be applied sparingly and excess material removed. Its barrier effect may impair the process of bone healing and obstruct the treatment of any concurrent infection.

## ADVERSE REACTIONS

Bone wax may cause a mild inflammatory reaction in the surrounding tissue and inhibit osteogenesis in the affected area.

## STERILITY

Bone wax is supplied sterile (sterilization method: gamma irradiation). Do not resterilize, critical changes in suitability for use possible (more details available on request)! Do not use damaged packages. Packages that have been opened for use but then not used are to be discarded.

## STORAGE CONDITIONS

Bone wax should be stored in a clean, dry place below 25° C. It should not be used after the expiry date.

## PRESENTATION

Bone wax is supplied sterile in packs of 24 x 2.5 gram slabs.

## EXPLANATION OF SYMBOLS ON THE PACKAGE

**CE** 1275 EC symbol and identification number of the notified body.  
The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

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