SERAFAST®

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**POLYGLYCOLIC ACID-CAPROLACTONE
SUTURE MATERIAL**
• SERAFAST®

DESCRIPTION

SERAFAST® is a synthetic absorbable monofilament sterile surgical suture material. It is composed of a segmented copolymer with the basic compounds glycolide and ε-caprolactone, empirical formula $(C_7H_{12}O_2)_m(C_6H_{10}O_2)_n$. It is available dyed with D+C violett Nr. 2, Colour Index Nr. 60725, or undyed.

The suture material complies with the underlying requirements of the Medical Devices Directive 93/42/EEC and the European Pharmacopoeia (EP) monograph "Sutures, sterile synthetic absorbable (Fila resorbilia synthetica monofilamenta sterilia)".

INDICATIONS

SERAFAST® is used to approximate soft tissues or for tying ligatures in cases where rapidly absorbable (up to 2 weeks) suture is indicated.

USE

When selecting and using the suture the condition of the patient, the surgeon's experience, the surgical procedure, and the type and size of the wound should be taken into account. Knot security requires the established standard knotting techniques for synthetic monofilament sutures. To avoid damage to the needle while suturing it is recommended that the needle is grasped towards the needle base in an area between the end of the first third and the middle of the needle.

ACTION

The purpose of the surgical suture is to create tissue-to-tissue approximation during the treatment of wounds. The ligature, or tie, as a special form of suture serves to occlude lumens. SERAFAST® is broken down in the tissue by hydrolysis to adipic acid which is then metabolised in the body. This absorption expresses itself as a loss

of tensile strength which is later accompanied by loss of mass. The tensile strength is reduced to 50% within approx. 8 – 13 days, and to 0% of the original breaking load after approx. 4 weeks. The absorption process is completed after about 90 – 120 days.

CONTRAINDICATIONS

SERAFAST® must not be used if longer-term or permanent stability of the suture material is required, e.g. securing synthetic implants (heart valves etc.).

WARNINGS / PRECAUTIONS / INTERACTIONS

The demands placed on the suture material and the required techniques will vary, depending on the area of application. The user should therefore be familiar with the various surgical techniques and procedures. Its in-vivo properties should be considered when selecting the suture material. The use of SERAFAST® should be carefully considered in patients with delayed wound healing. As with any foreign body there is a risk of calculus formation with prolonged contact with salt solutions (e.g. in the biliary or urinary tracts). Sutures under tension or those needing additional support may require the supplemental use of non-absorbable sutures or suitable immobilisation procedures.

In intracutaneous sutures the suture material – in particular the knots – should be placed as deep as possible. The surgical instruments should be handled as carefully as possible to avoid damage to the suture. Bending the needles or grasping them outside the recommended area may result in needle breakage. The handling and disposal of needles should be carried out with particular care due to the risk of contamination following needle stick injuries.

ADVERSE REACTIONS

Even with proper use of SERAFAST® the following adverse effects may occur: Transient localised irritation, transient inflammatory foreign-body reactions, calculus formation with prolonged contact with salt solutions.

STERILITY

SERAFAST® is supplied sterile (sterilization procedure: ethylene oxide gas). Do not re-sterilize, critical changes in suitability for use possible (more details available on request)! Do not use damaged packages. Discard opened unused packs!

Do not use if package is damaged

Do not re-use

STORAGE CONDITIONS

No special storage conditions are required. Do not use after the expiry date.

Sterilized with ethylene oxide

DISPOSAL INFORMATION

Dispose of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

PRESENTATION

SERAFAST® is available in various strengths and lengths, dyed and undyed, attached to various needle types, and unneeded. The sutures and their combinations are supplied sterile in single and multiple packs (multipack). Please see catalogue for precise details.

 1014 CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

Surgical suture, absorbable, monofilament, undyed

Surgical suture, absorbable, monofilament, dyed

Polyglycolic acid – ϵ -caprolactone

Antireflective needle

Detachable needle

Needle identification symbols and abbreviations are explained in the catalogue.

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SERAG-WIESSNER GmbH & Co. KG