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**POLYGLYCOLIC ACID-CAPROLACTONE
SUTURE MATERIAL**
• SERAFAST®

WARNINGS AND PRECAUTIONS

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified and trained personnel. It is the physician's responsibility to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: ethylene oxide).
- The product is intended for single use only.
- Do not resterilize or reuse the product as critical changes in performance and safety are possible.
- Do not use after expiration date.
- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or accidentally opened prior to use.
- Use the product only for the intended indication.
- Please contact the manufacturer for further information.

of needles made of surgical steel. SERAFAST® may be supplied in combination with other components, such as a suture guide tube. The suture material complies with the underlying requirements of the European Pharmacopoeia (EP) monograph "Sutures, sterile synthetic absorbable monofilament (Fila resorbilia synthetica monofilamenta sterilia)".

SERAFAST® is used to approximate soft tissues or for tying ligatures in cases where rapidly absorbable (up to 2 weeks) suture is indicated.

INSTRUCTIONS FOR USE

When selecting and using the suture the condition of the patient, the surgeon's experience, the surgical procedure, and the type and size of the wound should be taken into account. Knot security requires the established standard knotting techniques for synthetic monofilament sutures. To avoid damage to the needle while suturing it is recommended that the needle is grasped towards the needle base in an area between the end of the first third and the middle of the needle.

ACTION

The purpose of the surgical suture is to create tissue-to-tissue approximation during the treatment of wounds. The ligature, or tie, as a special form of suture serves to occlude lumens. SERAFAST® is broken down in the tissue by hydrolysis to adipic acid which is then metabolised in the body. This absorption expresses itself as a loss of tensile strength which is later accompanied by loss of mass. The breaking strength is reduced to about 60% of its initial value within 7 days (in-vitro). Absorption is completed after about 90 days (in-vivo).

CONTRAINDICATIONS

SERAFAST® must not be used if longer-term or permanent stability of the suture material is required, e.g. securing synthetic implants (heart valves etc.). The device must not be used in patients with allergies to any of the product components.

DESCRIPTION AND INTENDED PURPOSE

SERAFAST® is a synthetic absorbable monofilament surgical suture material. It is composed of a segmented copolymer with the basic compounds glycolide and ε-caprolactone. It is available dyed with D+C violet Nr. 2, Colour Index Nr. 60725, or undyed. SERAFAST® is available in various gauge sizes and lengths. It is supplied with various types

RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS

The demands placed on the suture material and the required techniques will vary, depending on the area of application. The user should therefore be familiar with the various surgical techniques and procedures. Its properties should be considered when selecting the suture material. The use of SERAFAST® should be carefully considered in patients with delayed wound healing. As with any foreign body there is a risk of calculus formation with prolonged contact with salt solutions (e.g. in the biliary or urinary tracts). Sutures under tension or those needing additional support may require the supplemental use of non-absorbable sutures or suitable immobilisation procedures.

In intracutaneous sutures the suture material – in particular the knots – should be placed as deep as possible. The surgical instruments should be handled as carefully as possible to avoid damage to the suture. Bending the needles or grasping them outside the recommended area may result in needle breakage. The handling and disposal of needles should be carried out with particular care due to the risk of contamination following needle stick injuries.

Considerations should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur. The use of absorbable sutures may be inappropriate in elderly, malnourished or debilitated patients. Skin sutures which remain in place longer than 7 days may cause local irritation and should be snipped off or removed.

ADVERSE REACTIONS

Even with proper use of SERAFAST® the following adverse effects may occur:

Transient localised irritation, transient inflammatory foreign-body reactions, calculus formation with prolonged contact with salt solutions.

Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

STORAGE CONDITIONS

The storage conditions apply as indicated on the label.

DISPOSAL INFORMATION

Disposal of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

MEANING OF THE SYMBOLS

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG