# SERAMON®

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 GB

## INSTRUCTIONS FOR USE

### SUTURE MATERIAL of POLYTETRAFLUOROETHYLENE

- SERAMON® undyed
- SERAMON® CHORDAE LOOPS

## DESCRIPTION

SERAMON® is a synthetically manufactured, monofilament, non-absorbable, sterile surgical suture material. The basic suture material consists of polytetrafluoroethylene (PTFE), chemical formula  $(C_2F_4)_n$ , which has been highly compressed and stretched to achieve a high tensile strength. SERAMON® is only available in undyed form. The PTFE suture materials are available as combination sutures with or without stainless steel needles and, if required, with specific accessories such as pledgets for cardiac and vascular sutures. They satisfy the essential requirements of the Medical Device Directive 93/42/EEC and apart from possible minimal deviations in the thread diameter, the harmonised standard of the European Pharmacopoeia "Sutures, sterile non-absorbable (Fila non resorbilia sterilia)".

## INDICATIONS

SERAMON® is indicated for soft-tissue adaptation or ligation, in which non-absorbable suture material is indicated, including use in the central circulatory and nervous systems. There has been insufficient experience with use in microsurgery

on peripheral vessels or in ophthalmology. PTFE threads can also be used as retaining sutures and for marking purposes.

## USE

When selecting and using suture material, the patient's condition, the physician's experience, the surgical technique used and the wound size should be taken into consideration. The established standard knot-tying techniques should be used to ensure knot-holding security. As a result of the even thread surface, additional knots may be necessary in tissues that are under tension. In order to prevent needle damage during suturing, it is recommended that the needle be grasped at a point between the end of its first third (as seen from its end) and the midpoint of the needle.

## ACTION

In wound care, the function of a medical suture is to attach tissue to tissue or tissue to allogeneic and/or xenogeneic material. Ligation or ligature as a special suturing form is used in the closure of hollow organs.

After an initially minimal tissue reaction, suture material made from PTFE is slowly and progressively encapsulated by connective tissue in the body. It is non-absorbable, physiologically inert and thus completely retains its tensile strength. As a result of its monofilament thread structure, SERAMON® can also be used in tissue at risk of infection.

## CONTRAINDICATIONS

not known.

## WARNINGS / PRECAUTIONS / INTERACTIONS

The demands on the suture material and the techniques required vary according to the therapeutic indication. The user should therefore be familiar with the surgical techniques. When selecting the suture material, its in-vivo properties should be taken into account. As with all foreign bodies, SERAMON® is associated with the risk of calculus formation following long-term contact with the bile ducts or urinary tract. For sutures placed in infected tissues, use of appropriate concomitant

measures is recommended. The surgical instruments should be handled carefully, in order to avoid damage to threads. Deformation of needles or holding the needle outside the recommended area can result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled and disposed of with particular care.

## ADVERSE REACTIONS

With the correct use of SERAMON®, the following adverse reactions can occur: transient local irritation, inflammatory reactions to the foreign body and, in rare cases, formation of suture granulomas or fistulas.

## STERILITY

SERAMON® is delivered sterile (sterilization method: ethylene oxide). Do not resterilize, critical changes in suitability for use possible (more details available on request)! Do not use damaged packs!

Packs that have been opened for use but then not used are to be discarded!

## STORAGE CONDITIONS

The suture material should be stored in a clean and dry place at room temperature. Do not use after the expiry date!

## PRESENTATION

SERAMON® is supplied in a range of gauge sizes and lengths in the form of needle-thread combinations with various needle types or without needles. The combination sutures are available in a sterile pack, individually or with several items in a multipack, as well as furnished with PTFE pledgets to reinforce the body's own structures. More precise details are given in the catalogue.

## EXPLANATION OF SYMBOLS ON THE PACKAGE

**CE**<sub>1275</sub> CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

Surgical suture, nonabsorbable, monofilament, undyed

Polytetrafluoroethylene

Antireflective needle

Detachable needle

Pledget (length times width times thickness [mm])

Needle identification symbols and abbreviations are explained in the catalogue.

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SERAG-WIESSNER GmbH & Co. KG