CE 1014

ce240009 / N01599v4

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SUTURE MATERIAL of POLYTETRAFLUOROETHYLENE

- **SERAMON® undyed**
- **SERAMON® CHORDAE LOOPS**

WARNINGS AND PRECAUTIONS

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified and trained personnel. It is the physician's responsibility to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: ethylene oxide). 
- The product is intended for single use only. 
- Do not resterilize or re-use the product as critical changes in performance and safety are possible. 
- Do not use after expiration date. 
- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or unintentionally opened prior to use. 
- Use the product only for the intended indication.

DESCRIPTION AND INTENDED PURPOSE

SERAMON® is a synthetic, monofilament, non-absorbable surgical suture material. The raw material for the suture is composed of polytetrafluoroethylene (PTFE). SERAMON® is supplied only in undyed form. It is available in various gauge sizes and lengths. The PTFE suture materials are supplied in combination with various types of needles made of surgical steel or without needles. Pledgets to reinforce the body's own structures may be enclosed with SERAMON®. Apart from marginally deviations with respect to the diameter, the suture material complies with the current edition of the European Pharmacopoeia monograph (EP) for "Sutures, sterile non-absorbable (Fila non resorbilia sterilia)".

SERAMON® is indicated for soft-tissue adaptation or ligation, in which non-absorbable suture material is indicated, including use in the central circulatory and nervous systems. PTFE threads can also be used as holding sutures and for marking purposes.

INSTRUCTIONS FOR USE

In the choice and use of suture material, due account should be taken of the condition of the patient, the experience of the physician, the surgical technique, and the type and size of the wound. Standard knot-tying techniques should be used to ensure knot-holding security. As a result of the even thread surface, additional knots may be necessary in tissues that are under tension. In order to prevent needle damage during suturing it is recommended that the needle be grasped at an area between the end of its first third (as seen from its end) and its midpoint.

Unneeded suture material may be used for ligation or may be connected with surgical eye needles. For a list of suitable needles and application instructions, please refer to the Instructions for Use of Serag-Wiessner sterile eye needles. The suitability of other

manu-facturers' products has not been tested. When using other brands, their specific instructions must be observed.

ACTION

The function of a medical suture in wound care is to attach tissue to tissue or tissue to foreign material. Ligatures are a special type of suture intended to close hollow organs. After an initially minimal tissue reaction, suture material made from PTFE is slowly and progressively encapsulated by connective tissue in the body. It is non-absorbable, physiologically inert and thus completely retains its tensile strength.

CONTRAINDICATIONS

The device must not be used in patients with allergies to any of the product components.

RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS

There has been insufficient experience with the use in microsurgery, ophthalmology and surgical interventions at peripheral nerves. The requirements of the suture material and the techniques needed vary with the indication. In the choice of suture material, account should be taken of the characteristics of the suture material concerned. Like all foreign bodies, the device can give rise to calculus formation when in prolonged contact with salt-containing solutions (e.g. in the biliary or urinary tract). As a result of its monofilament thread structure, SERAMON[®] can also be used in tissue at risk of infection. In the case of sutures in infected tissue the use of suitable ancillary measures is called for. Surgical instruments should be handled with care in order to avoid damage to thread. Bending of needles or grasping them outside of the recommended area can result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled with particular care.

ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of the device: local irritation, inflammatory reactions to the foreign body; formation of foreign-body giant cells, suture fistulas, or granulomas.

Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

STORAGE CONDITIONS

The storage conditions apply as indicated on the label.

DISPOSAL INFORMATION

Disposal of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

MEANING OF THE SYMBOLS

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG