PTFE PLEDGETS

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PTFE PLEDGETS

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PLEDGETY ZPTFE

en INSTRUCTIONS FOR USE

PTFE PLEDGETS

tr KULLANIM TALİMATI

PTFE PLEJİTLER

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PLEDGET IN PTFE

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PLEDGETS DE PTFE

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PTFE PLEDGETS

DESCRIPTION

PTFE Pledgets (suture buttresses) are made of non-absorbable, undyed polytetrafluorethylene (PTFE) polymer fibres, which have been compressed to form a felt. The PTFE Pledget thus has a stable structure with a minimal pore size, small volume and high suture pullout strength.

INDICATIONS

The main function of the pledgets is to provide durable tissue reinforcement and to form a suture buttress in particularly vulnerable regions, especially the heart and vessels. By providing reinforcement for tissues or prostheses, the pledgets render suture pullout difficult from tissues subject to strong mechanical loading.

USE

The pledgets are designed for use in combination with surgical suture material. They generally have prepunched holes to enable the user to thread them easily. The patient's condition, the experience of the surgeon and the surgical technique should be considered when selecting and using the pledgets.

ACTION

PTFE caused only minimal tissue reactions during clinical studies. PTFE Pledgets are neither absorbed nor weakened by tissue enzymes; nor are they degraded in the presence of an infection.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS / PRECAUTIONS / INTERACTIONS

Improper handling with surgical instruments or contact with sharp edges should be avoided to prevent fraying of fibres or other mechanical damage. PTFE Pledgets should not be exposed to laser beams. Perforating the pledgets outside the

prepunched holes may reduce needle insertion force.

ADVERSE REACTIONS

Even with proper use of PTFE Pledget, the following undesired effects may occur: transient local irritation, inflammatory foreign-body reactions, fistula or granuloma formation.

STERILITY

PTFE Pledgets have been sterilised using ethylene oxide gas. Do not resterilize, critical changes in suitability for use possible (more details available on request)! Do not use damaged packs! Discard unused opened packs!

STORAGE CONDITIONS

PTFE Pledgets should be stored in a clean dry place at room temperature. Do not use after the expiry date!

PRESENTATION

Please refer to the catalogue for exact details.

EXPLANATION OF SYMBOLS ON THE PACKAGE

CE₁₀₁₄ CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

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SERAG-WIESSNER GmbH & Co. KG