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- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified and trained personnel. It is the physician's responsibility to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: ethylene oxide).
- The product is intended for single use only.
- Do not resterilize or reuse the product as critical changes in performance and safety are possible.
- Do not use after expiration date.
- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or accidentally opened prior to use.
- Use the product only for the intended indication.
- Please contact the manufacturer for further information.

STERILE EO

**DESCRIPTION AND INTENDED PURPOSE**

SERAFIT® PROTECT is a synthetic, absorbable, multifilament braided surgical suture material. It is made from the homopolymer of glycolic acid. SERAFIT® PROTECT is coated with calcium stearate, polyglycolic acid/poly lactide and chlorhexidine diacetate. It is supplied dyed with D+C violet no. 2, Colour Index No. 60725 or undyed. SERAFIT® PROTECT is available in various gauge

sizes and lengths. It is supplied with various types of needles made of surgical steel or without needles.

The suture material complies with the European Pharmacopoeia (EP) monograph "Sutures, sterile synthetic absorbable braided (Fila resorbilia synthetica torta sterilia)".

SERAFIT® PROTECT is intended for use in adaptation of soft tissues and in ligatures where medium-term absorbable thread material is required.

INSTRUCTIONS FOR USE

In the choice and use of suture material, due account should be taken of the condition of the patient, the experience of the doctor, the surgical technique, and the type and size of the wound. Standard knot-tying techniques should be used to ensure knot-holding security. In order to prevent needle damage during suturing it is recommended that the needle be grasped at an area between the end of its first third (as seen from its end) and its midpoint.

Unneeded suture material may be used for ligation or may be connected with surgical eye needles. For suitable needles and application instructions, please refer to the Instructions for Use of Serag-Wiessner's sterile eye needles. The suitability of other manufacturers' products has not been tested. When using other brands, their specific instructions must be observed.

ACTION

The function of a medical suture in wound care is to attach tissue to tissue. Ligatures are a special type of suture intended to close hollow organs. SERAFIT® PROTECT is degraded by hydrolysis in tissues to form glycolic acid, which is then metabolized in the body. This absorption becomes manifest firstly as a reduction in tensile strength and later as a loss of mass of the thread. The breaking strength falls to about 60% of its initial value after 14 days (in-vitro). Absorption is completed after about 90 days (in-vivo). In vitro and in vivo studies with SERAFIT® PROTECT have demonstrated that colonisation of the suture mate-

rial with *Staphylococcus aureus*, *Staphylococcus epidermidis* incl. the MRSA and MRSE variants is prevented or reduced. The clinical efficacy of this property has not been demonstrated.

CONTRAINDICATION

SERAFIT® PROTECT must not be used where more prolonged or long-term stability of the suture material is required. Do not use SERAFIT® PROTECT in patients with known allergy or hypersensitivity to chlorhexidine or other components of the device. It is not intended to be used in direct contact with the heart, the central circulatory system or the central nervous system. Due to inadequate data, the application in microsurgery, orthopaedics, ophthalmology and otology or to other sensitive organs and newborns is not recommended.

RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS

The requirements of the suture material and the techniques required vary with the indication. Therefore, the user should be familiar with the surgical techniques to be used. In the choice of suture material, account should be taken of the characteristics of the suture material concerned. In patients with delayed wound healing SERAFIT® PROTECT should be used only after careful consideration. Considerations should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur. The use of absorbable sutures may be inappropriate in elderly, malnourished or debilitated patients. Like all foreign bodies, SERAFIT® PROTECT can give rise to calculus formation when in prolonged contact with salt-containing solutions (e.g. in the biliary or urinary tract). Additional use of nonabsorbable suture materials or suitable immobilization techniques may be required with sutures that are under tension or require greater support.

Use of SERAFIT® PROTECT does not preclude the need for compliance with generally applicable hygiene guidelines. Use of suitable ancillary measures is required with sutures in infected tissue. Because of the increased risk of infection as-

sociated with the multifilament structure of the thread, skin sutures should if possible be intracutaneous and the suture material, in particular the knots, should be placed as deep as possible. Skin sutures which remain in place longer than 7 days may cause local irritation and should be snipped off or removed.

Surgical instruments should be handled with care in order to avoid damage to the thread. Bending of needles or grasping them outside of the recommended area can result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled and disposed of with particular care.

ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of SERAFIT® PROTECT: transient local irritation, transient inflammatory reactions to the foreign body; rarely, formation of suture fistulas or granulomas, formation of calculus after prolonged contact with salt-containing solutions.

Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

STORAGE CONDITIONS

The storage conditions apply as indicated on the label.

DISPOSAL INFORMATION

Disposal of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

MEANING OF THE SYMBOLS

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG