**SERACOR® SULENE® TERYLENE**

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- **SERACOR<sup>®</sup>**
- **SULENE<sup>®</sup>**
- **Terylene**

### DESCRIPTION

SERACOR<sup>®</sup>, SULENE<sup>®</sup>, and Terylene are synthetic, braided, nonabsorbable, sterile surgical suture materials of polyethylene terephthalate (PET). PET is a polymer with the molecular formula  $(C_{10}H_8O_4)_n$ , which is generally synthesized from a terephthalic acid derivative and ethylene glycol. The suture materials are supplied in combinations with stainless steel needles or without needles. They comply with the underlying requirements of the Medical Devices Directive 93/42/EEC and the European Pharmacopoeia monograph "Sutures, sterile non-absorbable" (fila non resorbilia sterilis) and - concerning gauge size- the monograph "Nonabsorbable surgical suture" of the United States Pharmacopeia (USP) respectively.

- **SERACOR<sup>®</sup>** green and undyed  
SERACOR<sup>®</sup> is a precision-braided PET suture with a special silicone coating. SERACOR<sup>®</sup> green is dyed with D & C No. 6, colour index no. 61565.
- **SULENE<sup>®</sup>** green  
SULENE<sup>®</sup> is a precision-braided PET suture with a special silicone coating. It is dyed with D & C No. 6, colour index no. 61565.
- **Terylene** green and undyed  
Terylene is a braided, silicone-coated PET suture. Terylene green is dyed with D & C No. 6, colour index no. 61565.

### INDICATIONS

SERACOR<sup>®</sup>, SULENE<sup>®</sup>, and Terylene are suitable for adaptation of soft tissues and for ligatures where nonabsorbable suture material is indicated, including use in the central vascular and nervous systems and in particular where long-term stability of the suture material is required. Because of its special braiding and coating, SERACOR<sup>®</sup> is especially suitable for use in the central vascular

system. Polyester threads can also be used as holding sutures and for marking purposes.

### USE

In the choice and use of suture material, due account should be taken of the condition of the patient, the experience of the doctor, the surgical technique, and the size of the wound. Standard knot-tying techniques should be used to ensure knot-holding security. In order to prevent needle damage during suturing it is recommended that the needle be grasped at a point between the end of its first third (as seen from its end) and its midpoint.

### ACTION

The function of a medical suture in wound care is to attach tissue to tissue or tissue to allogeneic material. Ligatures are a special type of suture used to close hollow organs.

In the body, PET suture material is slowly and progressively encapsulated by connective tissue. It is nonabsorbable and does not lose tensile strength in tissues even after prolonged periods.

### CONTRAINDICATIONS

None known.

### WARNINGS / PRECAUTIONS / INTERACTIONS

The requirements of the suture material and the techniques needed vary with the indication. Therefore, the user should be familiar with the surgical techniques to be employed. In the choice of suture material, account should be taken of the in-vivo characteristics of the suture material concerned. Like all foreign bodies, polyester threads can give rise to calculi when in prolonged contact with the biliary or urinary tract. In the case of sutures in infected tissue the use of suitable ancillary measures is called for, as the multifilament thread structure can prolong infection. Surgical instruments should be handled with care in order to avoid damage to threads. Bending of needles or grasping them outside of the recommended area can result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled and disposed of with particular care.

## ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of polyester threads: local irritation and inflammatory reactions to the foreign body; formation of foreign-body giant cells, suture fistulas, or granulomas.

## STERILITY

SERACOR®, SULENE®, and Terylene are supplied sterile (sterilization method: ethylene oxide). Do not resterilize, critical changes in suitability for use possible (more details available on request)! Do not use damaged packs! Packs that have been opened for use but then not used are to be discarded!

## STORAGE CONDITIONS

No special storage conditions are required. Do not use after the expiry date.

## INSTRUCTIONS FOR REMOVAL FROM THE PACK

**Minibox:** For easy removal the thread should be withdrawn not obliquely, but vertically upwards out of the minibox.

## DISPOSAL INFORMATION

Dispose of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

## PRESENTATION

Our PET suture materials are supplied in various gauge sizes, lengths, and colours and with various types of needle or without needles. Pledgets to reinforce the body's own structures may be enclosed with SERACOR®. The threads and combinations are supplied in individual packs and multiple packs in a sterile package. More precise details are given in the catalogue.

## EXPLANATION OF SYMBOLS ON THE PACKAGE

 1014 CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

Surgical suture, nonabsorbable, braided, coated, undyed

Surgical suture, nonabsorbable, braided, coated, dyed

Polyester

Antireflective needle

Detachable needle

Pledget (length times width times thickness [mm])

Needle identification symbols and abbreviations are explained in the catalogue.

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SERAG-WIESSNER GmbH & Co. KG