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- SERACOR®
- SULENE®
- Terylene

WARNINGS AND PRECAUTIONS

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified and trained personnel. It is the physician's responsibility to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: ethylene oxide).
- The product is intended for single use only.
- Do not resterilize or reuse the product as critical changes in performance and safety are possible.
- Do not use after expiration date.
- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or accidentally opened prior to use.
- Use the product only for the intended indication.
- Please contact the manufacturer for further information.

colour index no. 61565 or undyed. The PET suture materials are available in various gauge sizes and lengths. They are supplied in combinations with various types of needles made of surgical steel or without needles. Pledgets may be enclosed with SERACOR® to reinforce the body's own structures. SERACOR®, SULENE® und Terylene comply with the current edition of the European Pharmacopoeia monograph "Sutures, sterile non-absorbable" (fila non resorbilia sterilia)".

SERACOR®, SULENE®, and Terylene are suitable for adaptation of soft tissues and for ligatures where nonabsorbable suture material is indicated, including use in the central vascular and nervous systems and in particular where long-term stability of the suture material is required. Polyester threads can also be used as holding sutures and for marking purposes.

INSTRUCTIONS FOR USE

In the choice and use of suture material, due account should be taken of the condition of the patient, the experience of the doctor, the surgical technique, and the size of the wound. Standard knot-tying techniques should be used to ensure knot-holding security. In order to prevent needle damage during suturing it is recommended that the needle be grasped at a point between the end of its first third (as seen from its end) and its midpoint.

Unneedled suture material may be used for ligation or may be connected with surgical eye needles. For suitable needles and application instructions, please refer to the Instructions for Use of Serag-Wiessner sterile eye needles. The suitability of other manufacturers' products has not been tested. When using other brands, their specific instructions must be observed.

ACTION

The function of a medical suture in wound care is to attach tissue to tissue or tissue to allogeneic material. Ligatures are a special type of suture used to close hollow organs.

In the body, PET suture material is slowly and pro-

DESCRIPTION AND INTENDED PURPOSE

SERACOR®, SULENE®, and Terylene are synthetic, braided, nonabsorbable surgical suture materials made of polyethylene terephthalate (PET).

The suture materials are coated with silicone. They are supplied dyed green with D&C No. 6,

gressively encapsulated by connective tissue. It is nonabsorbable and does not lose tensile strength in tissues even after prolonged periods.

CONTRAINDICATIONS

The device must not be used in patients with allergies to any of the product components.

RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS

The requirements of the suture material and the techniques needed vary with the indication. Therefore, the user should be familiar with the surgical techniques to be employed. In the choice of suture material, account should be taken of the characteristics of the suture material concerned. Like all foreign bodies, polyester threads can give rise to calculi when in prolonged contact with the biliary or urinary tract. In the case of sutures in infected tissue the use of suitable ancillary measures is called for, as the multifilament thread structure can prolong infection. Surgical instruments should be handled with care in order to avoid damage to threads. Bending of needles or grasping them outside of the recommended area can result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled and disposed of with particular care.

ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of polyester threads: local irritation and inflammatory reactions to the foreign body; formation of foreign-body giant cells, suture fistulas, or granulomas.

Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

STORAGE CONDITIONS

The storage conditions apply as indicated on the label.

INSTRUCTIONS FOR REMOVAL FROM THE PACK

Minibox: For easy removal the thread should be

withdrawn not obliquely, but vertically upwards out of the minibox.

DISPOSAL INFORMATION

Dispose of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

MEANING OF THE SYMBOLS

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG