POLYESTERBAND

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en

INSTRUCTIONS FOR USE

**POLYESTER TAPE,
braided, nonabsorbable,
undyed or dyed, sterile**

- Shirodkar cervical suture
- Polyester tape (unneeded)

DESCRIPTION

The polyester tape is a synthetic, braided, nonabsorbable, sterile tape composed of polyethylene terephthalate (PET). PET is a polymer with the empirical formula $[C_{10}H_8O_4]$, and is normally synthesized from a terephthal acid derivative and ethylene glycol. It is available dyed with D+C Green, Colour Index Nr. 61565, or undyed. Polyester tapes are available with needles for use as Shirodkar cervical sutures or unneeded.

INDICATIONS

Shirodkar cervical sutures are intended for use in the Shirodkar cervical cerclage procedure. The unneeded polyester tapes are used as stay or retracting tapes.

USE

Armed with a curved needle the polyester tape is inserted beneath the cervical mucosa until it encircles the entire cervix. The tape is then knotted and tightened to close the cervix, yet leaving a slight opening. In order to prevent damage to the needle while suturing it is recommended that the needle be grasped towards the swaged end at a point between the end of the first third and the mid-point.

The unneeded polyester tapes can be used to temporarily retract organs and parts of organs out of the surgical field, thus facilitating the surgeon's work in that area. Afterwards the tapes are removed.

ACTION

PET tapes are nonabsorbable and do not lose their tensile strength in the tissue, even over prolonged periods of time. The Shirodkar cervical

suture adequately reinforces closure of the cervix in the presence of premature opening or cervical incompetence during the second trimester of pregnancy. The cerclage is removed at the end of the pregnancy.

The unneeded polyester tapes are suitable for use as slings to retract organs and other structures because their flat form and pliable nature prevents them from cutting into the tissue.

CONTRAINDICATIONS

None known.

WARNINGS / PRECAUTIONS / INTERACTIONS

At the end of pregnancy the cervical cerclage should be removed to avoid erosions into neighbouring tissues and organs. The surgical instruments used should be handled with care to avoid damage to the polyester tape. Bending the needles or grasping them outside the recommended area can result in needle breakage. Needles should be handled and disposed of with particular care due to the risk of infection from needle stick injuries.

ADVERSE REACTIONS

The following adverse reactions may occur despite correct use of the polyester tape:

Minimal transient tissue reaction in the form of localized irritation and inflammatory foreign-body reactions; infection may be prolonged when used in infected tissues.

STERILITY

The polyester tapes are gas sterilised with ethylene oxide. Do not resterilize, critical changes in suitability for use possible (more details available on request)! Do not use damaged packs! Packs that have been opened for use but then not used are to be discarded!

STORAGE CONDITIONS

The packets should be stored in a clean dry place at room temperature. Do not use after the expiry date!

PRESENTATION

The Shirodkar cervical suture is available together with a variety of needles and in various lengths. The unneeded polyester tapes are also available in various lengths. All products are supplied sterile in single packets. Please see catalogue for precise details.

EXPLANATION OF SYMBOLS ON THE PACKAGE

CE₁₂₇₅ CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

Needle identification symbols and abbreviations are explained in the catalogue.

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SERAG-WIESSNER GmbH & Co. KG