SERAFLEX®

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SILK SUTURE MATERIALS**• SERAFLEX®****DESCRIPTION**

SERAFLEX® is a natural, nonabsorbable, sterile surgical suture material derived from the cocoon fibres of the silkworm moth (*Bombyx mori* L.). The starting material, the raw silk fibre, consists mostly of the proteins fibroin and sericin.

In the manufacture of SERAFLEX® the raw silk fibre is degummed, i.e. freed of sericin, and braided. SERAFLEX® is a braided, silicone-coated suture and is available dyed black with Logwood black (haematin), Colour Index No. 75290 or undyed. SERAFLEX® in the cassette pack is supplied in an isopropanol-containing preserving solution.

Our silk suture materials are supplied in combinations with stainless steel needles or without needles. They comply with the underlying requirements of the Medical Device Directive 93/42/EEC and the current edition of the European Pharmacopoeia monograph (EP) for "sutures, sterile non-absorbable (fila non resorbilia sterilia)" or of the monograph "Nonabsorbable surgical suture" of the United States Pharmacopoeia (USP). In deviation from this standard, the maximum thread length in the case of the cassette packs is not 4 but 100 m.

INDICATIONS

Our silk suture materials are suitable for all wound sutures and ligatures in which nonabsorbable thread material is indicated with the exception of use in the central vascular or nervous systems. Silk threads can also be used as holding sutures and for marking purposes.

USE

In the choice and use of suture material, due account should be taken of the condition of the patient, the experience of the doctor, the surgical

technique, and the size of the wound. Standard knot-tying techniques should be used to ensure knot-holding security. In order to prevent needle damage during suturing it is recommended that the needle be grasped at a point between the end of its first third (as seen from its end) and its midpoint.

ACTION

The function of a medical suture in wound care is to attach tissue to tissue or tissue to allogeneic material. Ligatures are a special type of suture intended to close hollow organs.

In the body, silk suture material is slowly and progressively encapsulated by connective tissue. Although it is nonabsorbable, it can gradually lose tensile strength in tissues as a result of proteolytic breakdown.

CONTRAINDICATIONS

Our silk suture materials should not be used in persons known to be allergic to silk. They are not intended for use in the central vascular or nervous systems. Because of the loss of tensile strength that may occur in the body, silk suture materials should not be used when prolonged stability of the thread is required.

WARNINGS / PRECAUTIONS / INTERACTIONS

The requirements of the suture material and the techniques needed vary with the indication. Therefore, the user should be familiar with the surgical techniques to be employed. In the choice of suture material, account should be taken of the in-vivo characteristics of the suture material concerned. Like all foreign bodies, silk threads can give rise to calculus formation when in prolonged contact with the biliary or urinary tract. In the case of sutures in infected tissue the use of suitable ancillary measures is called for, as the multifilament thread structure can prolong infection. Surgical instruments should be handled with care in order to avoid damage to threads. Bending of needles or grasping them outside of the recommended area can result in needle breakage. Because of the risk of infection from needle

pricks, needles should be handled and disposed of with particular care.

ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of silk threads: moderately severe transient local irritation and inflammatory reactions to the foreign body; formation of suture fistulas or granulomas; formation of calculi after prolonged contact with salt-containing solutions; in predisposed patients, allergic reactions.

STERILITY

Our silk suture materials are gas-sterilized with ethylene oxide and Silk in the cassette pack is sterilized by irradiation. Do not resterilize, critical changes in suitability for use possible (more details available on request)! Do not use damaged packs!

For cassette packs an antiseptic removal after opening is guaranteed for a period of 6 months. Please note the expiry date on the outer pack by using the enclosed label.

Packs that have been opened for use but then not used — or in the case of the cassette packs not handled correctly (see INSTRUCTIONS FOR REMOVAL FROM THE PACK) — are to be discarded!

STORAGE CONDITIONS

No special storage conditions are required. Do not use after the expiry date or, with cassette packs, after the use by date! Store cassette packs protected from light and UV-radiation.

DISPOSAL INFORMATION

Dispose of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

INSTRUCTIONS FOR REMOVAL FROM THE PACK

Minibox: For easy removal the thread should be withdrawn not obliquely, but vertically upwards out of the minibox.

Minibox - microsutures in protective sleeve: Remove the sleeve from the minibox attachment;

grasp the needle with the needle holder and withdraw the thread from the protective sleeve; if a **short thread** is required, it should be cut through together with the protective sleeve; in the case of **double-armed threads**, two individual threads can, if required, be withdrawn after the loop of thread at the end of the protective sleeve has been cut through.

Cassette pack (also called suture dispenser or outer pack) : Remove the tamperproof closure (dispenser cap closure) and use the enclosed label to note the expiry date on the outer pack; observe the following instructions to guarantee aseptic removal for a period of 6 months after opening. Flip up the dispenser cap at the clip and fix it to the back part of the cap. If the outer pack has already been opened, i.e. the dispenser cap closure has been removed or opened, **disinfect** the inner surface of the dispenser **before use** with a commercially available alcohol-based product. Pull out the end of the suture with a sterile instrument. Discard the first 10 cm of the suture. Pull out the desired length of suture, cut it off flush, and close the cap. **After finishing the procedure, but no more than 4 hours later, or if required again disinfect** the inner surface of the dispenser and the cap with a commercially available alcohol-based product and close the cap.

PRESENTATION

Our silk suture materials are supplied in various gauge sizes, lengths, and colours and with various types of needle or without needles. The threads and combinations are supplied in individual packs and multiple packs in a sterile package. Economical cassette packs are available for bulk purchasers. More precise details are given in the catalogue.

EXPLANATION OF SYMBOLS ON THE PACKAGE

CE1014 CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

Surgical suture, nonabsorbable,
braided, coated, undyed

Surgical suture, nonabsorbable,
braided, coated, dyed

Silk

Antireflective needle

Detachable needle

Needle identification symbols and abbreviations
are explained in the catalogue.

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SERAG-WIESSNER GmbH & Co. KG