SERAMESH® SERATOM® SERATIM® SERASIS® SERATEX® SERAGYN®

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Sterile synthetic surgical mesh implant

DESCRIPTION

SERAG-WIESSNER textile mesh implants are knitted from synthetically produced polypropylene threads. These ground threads may be multifilament, monofilament or partially absorbable bi-component threads with an absorbable coating of glycolide and ϵ -caprolactone. Surgical suture materials may be attached to the mesh to act as stay sutures. The implants may be coloured with D+C violet no. 2, Colour Index no. 60725, Pthalocyanine blue, Colour Index no. 74160, Solvent blue, Colour Index no. 61568 or Pigment black 11, Colour Index no. 77499, over the whole surface of the mesh or in stripes for orientation.

The basic structure can be determined from the combination of letters given in the product description (see below, additional product labelling). Mesh implants made of non-absorbable polypropylene threads, which remain in the body, have shown unchanged stability properties in vivo over a long period of time. Scars form during the healing process, especially at the margins of the implant, and cells are able to grow into the pores and even completely penetrate and cover the textile implant.

The different knitted structures provide porous tear-resistant mesh implants of high tensile strength and with a flexibility that is appropriate to the anatomical areas for which they are intended. We can provide you with special data sheets on the individual designs on request.

INDICATIONS

SERAG-WIESSNER textile implants are intended to provide surgical support, act as a bridging material and/or provide reinforcement for the body's natural structures.

USE

Different dimensions and constructions are available for the various indications.

The textile implants are sold under the following brand names:

SERAMESH[®], SERASIS[®], SERATOM[®], SERATIM[®], SERATEX[®], SERAGYN[®]

ACTION

SERAG-WIESSNER textile implants reinforce or bridge damaged anatomical structures permanently.

CONTRAINDICATIONS

If mesh implants are used in younger patients who have not finished growing, or in women who may become or are planning to become pregnant, it must be remembered that mesh implants stretch to only a limited extent and do not grow with the patient.

The operation should not be performed during pregnancy. Otherwise, all the usual contraindications to (mesh implant) surgery apply.

WARNINGS / PRECAUTIONS / INTERACTIONS

Preoperative:

The surgeon should be familiar with the appropriate surgical techniques, especially those of minimally invasive surgery. Studies have shown that the main danger of erosion or other complications (such as bleeding, dyspareunia, renewed or increased incontinence) lies in the operative procedure. For this reason, we offer various workshops. Please contact us about them.

Intraoperative:

The implantation should be carried out with due care and attention in order to prevent any damage occurring to the major blood vessels, nerves and internal organs. Paying close attention to the local anatomy, as well as the correct passage of the instruments, minimises the risks.

Don't use a mesh implant in patients who are on anticoagulants.

If a mesh is implanted into infected tissue, the infection may be prolonged and there is a danger of fistula formation or rejection reactions. This applies particularly to the multifilament implants SERAMESH[®] SE and LE. It also has to be remembered that in the case of infection, the material may have to be removed.

Do not use urological or gynaecological pelvic floor implants in patients with urinary tract infections, in immunocompromised patients, in patients with renal failure or relative obstruction of the upper urinary tract.

The extremely light **SERAMESH**[®] XO should not be used to repair structures, such as large abdominal wall hernias, where considerable pressure will be exerted on the mesh.

It must be ensured that the mesh implant is inserted with the correct tension (usually tension-free, although **SERATIM**[®] is placed under tension), not twisted and wrinkle-free (especially while implantation of **SERAGYN**[®]).

Particular attention must be paid to the correct positioning, because if the mesh is put in the wrong position, it will be necessary to revise the implant or even remove it altogether and put a new one in the right place, as there is otherwise a risk of subsequent damage and late complications (e.g. mucosal erosion). The incontinence implants **SERASIS**[®] and **SERATIM**[®] allow adjustment of the tension of the tape until at least the second postoperative day with the aid of the tape ends or special stay sutures. If the implant is fitted individually, it must be checked whether the integrity of the implant is affected. Cystoscopy is advised in the case of operations near the bladder, in order to confirm that the bladder is still intact or to identify any perforation. This also applies to the intestine and other hollow organs. Adhesions may form if the implant comes into direct contact with the intestinal tract.

Postoperative:

Postoperative bleeding may occur. Check for any symptoms and signs of bleeding before discharging the patient from hospital.

After pelvic surgery, the patient should be told not to lift heavy objects or participate in any strenuous sports (e.g. cycling, jogging, swimming with breast-stroke leg movements) for at least 3-8 weeks (depending on the complexity of the operation) and to avoid sexual intercourse for a month. If dysuria, bleeding or any other problem should occur, the patient must contact the surgeon immediately.

Surgical instruments:

Even with normal use, all surgical instruments may show signs of wear or damage. Reusable instruments should be checked with the naked eye before use. Any instruments that have defects or are corroded must not be used, but are to be discarded.

The surgical instruments should be handled carefully in order to avoid damaging the mesh implant or any suture material that may have been used for fixation. The tapes, attachments or the stay sutures on the textile implants must be introduced into the eyelet of the instrument carefully so as to prevent individual thread systems from being threaded on, which could result in the tape or attachment tearing off. Textile implants may be mechanically damaged by clamps, staples or clips. The mesh implants must not be secured with clips or sutures at the margin – there must be a distance of at least 6 mm or 2 pores from the edge of the mesh.

ADVERSE REACTIONS

The following adverse reactions may occur even with the correct use of textile implants:

Mild transient tissue reactions in the form of local irritation and inflammatory foreign body reactions and the formation of seromas. Like all foreign bodies, textile implants can have a detrimental effect on pre-existing infections (see Warnings / Precautions / Interactions).

Exception for specifically labelled products, interactions of the textile implants in the MRI field cannot be observed yet. If the mesh implant is put under too much tension, e.g. after overcorrection, it may result in transient or permanent stenosis of the urinary tract or functional impairment of voiding the bladder or bowel.

Passage of the instruments may damage blood vessels, nerves, or hollow organs, causing perforations or tears that require surgical repair.

STERILITY

SERAG-WIESSNER textile implants are supplied in sterile packs (sterilisation method: ethylene oxide). Do not resterilize, critical changes in suit-

ability for use possible (more details available on request)! Do not use damaged packs! Discard unused opened packs!

STORAGE CONDITIONS

Store SERAG-WIESSNER textile implants in a clean dry place at room temperature. Do not use after the expiry date.

PRESENTATION

SERAG-WIESSNER textile implants are supplied in various designs and sizes. If required, they can be supplied with an incision to facilitate their adaptation to the anatomical structures. Sizes smaller than 15 x 15 cm are supplied flat, the rest are packed folded. For your patient's operation records, a product label with all data relevant to documentation is supplied with every mesh implant. An implant ID-card may be requested for the patient. Please consult our catalogue for more information.

EXPLANATION OF SYMBOLS ON THE PACKAGE


CE₁₂₇₅ CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

PP Polypropylene

PP/PGACL Polypropylene/
polyglycolic acid – ε-caprolactone

 Textile mesh implant with different structures, stitch widths and dimensions.

The symbol within the circle shows the basic structure of the mesh implant
In detail this can be:

 Monofilament, nonabsorbable, not dyed

 Monofilament, nonabsorbable, dyed

 Multifilament, nonabsorbable, not dyed

 Bicomponent fibre, partially absorbable, not dyed

 Bicomponent fibre, partially absorbable, dyed

ADDITIONAL IDENTIFICATION CODES:

SO stitch width < 2 mm, monofilament thread structure
SE stitch width < 2 mm, multifilament thread structure
LE stitch width > 4 mm, multifilament thread structure
XO stitch width > 2 mm, special monofilament structure
PA partially absorbable
MR visible in the MRI field

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SERAG-WIESSNER GmbH & Co. KG