# SULENE® GROSSPACKUNG TERYLENE

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- SULENE®
- TERYLENE

### DESCRIPTION

SULENE®, and TERYLENE are synthetic, braided, nonabsorbable, sterile surgical suture materials of polyethylene terephthalate (PET). PET is a polymer with the molecular formula  $(C_{10}H_8O_4)_n$ , which is generally synthesized from a terephthalic acid derivative and ethylene glycol.

The suture materials are supplied without needles. They satisfy the essential requirements of Directive 93/42/EEC on medical devices and the current edition of the harmonized standard of the European Pharmacopoeia for "sterile, nonabsorb-

able suture material (fila non resorbilia sterilia)" and - concerning gauge size- the monograph "Nonabsorbable surgical suture" of the United States Pharmacopeia (USP) respectively. In deviation from this standard, the maximum thread length is not 3.5 but 100 m.

- **SULENE® green**  
SULENE® is a precision-braided PET suture with a special silicone coating. It is dyed with D & C No. 6, colour index no. 61565.
- **TERYLENE green and undyed**  
TERYLENE is a braided, silicone-coated PET suture. TERYLENE green is dyed with D & C No. 6, colour index no. 61565.

### INDICATIONS

SULENE®, and TERYLENE are suitable for adaptation of soft tissues and for ligatures where non-absorbable suture material is indicated, in particular where long-term stability of the suture material is required. Polyester threads can also be used as holding sutures and for marking purposes.

### USE

In the choice and use of suture material, due account should be taken of the condition of the patient, the experience of the doctor, the surgical technique, and the size of the wound. Standard knot-tying techniques should be used to ensure knot-holding security.

### ACTION

The function of a medical suture in wound care is to attach tissue to tissue or tissue to allogeneic material. Ligatures are a special type of suture used to close hollow organs.

In the body, PET suture material is slowly and progressively encapsulated by connective tissue. It is nonabsorbable and does not lose tensile strength in tissues even after prolonged periods.

### CONTRAINDICATIONS

SULENE®, and TERYLENE in the cassette pack are not intended for use in the central vascular or nervous systems.

### WARNINGS / PRECAUTIONS / INTERACTIONS

The requirements of the suture material and the techniques needed vary with the indication. Therefore, the user should be familiar with the surgical techniques to be employed. In the choice of suture material, account should be taken of the in-vivo characteristics of the suture material concerned. Like all foreign bodies, polyester threads can give rise to calculi when in prolonged contact with the biliary or urinary tract. In the case of sutures in infected tissue the use of suitable ancillary measures is called for, as the multifilament thread structure can prolong infection. Surgical instruments should be handled with care in order to avoid damage to threads.

### ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of polyester threads: local irritation and inflammatory reactions to the foreign body; formation of foreign-body giant cells, suture fistulas, or granulomas.

### STERILITY

SULENE®, and TERYLENE are sterilized by irradiation. Do not resterilize, critical changes in suitability for use possible (more details available on request)! Do not use damaged packs!

After opening cassette packs aseptic removal is guaranteed for a period of 6 months. Please make a note of the use by date using the attached pack label.

Packs that have been opened for use but then not used or not handled correctly (see INSTRUCTIONS FOR REMOVAL FROM THE PACK) are to be discarded!

### STORAGE CONDITIONS

The suture material should be kept dry and clean at room temperature. Do not use after the expiry date or after the use by date respectively!

### INSTRUCTIONS FOR REMOVAL FROM THE PACK

Flip up the dispenser cap and pull out the end of the suture with a sterile instrument. Discard the first 10 cm of the suture. Pull out the desired

length of suture, cut it off flush, and close the cap. After finishing the procedure, but no more than 4 hours later, or if required again disinfect the inner surface of the dispenser and the cap with a commercially available alcohol-based product and close the cap.

## PRESENTATION

Our PET suture materials are supplied in various gauge sizes, lengths, and colours without needles. More precise details are given in the catalogue.

## EXPLANATION OF SYMBOLS ON THE PACKAGE

**CE** 1275 CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

Surgical suture, nonabsorbable, braided, coated, undyed

Surgical suture, nonabsorbable, braided, coated, dyed

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SERAG-WIESSNER GmbH & Co. KG