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POLYDIOXANONE SUTURE MATERIALS

- SERASYNTH®
- SERASYNTH® LOC
- SERASYNTH® SLIMSILING

DESCRIPTION

They are synthetically manufactured absorbable, monofilament, sterile surgical suture materials. They are manufactured of the polymer poly-p-dioxanone (PDO) with the formula $[C_6H_8O_3]_n$. It is available dyed with D&C violet Nr. 2, Colour Index Nr. 60725. An anchor fixation system based on monidirectional grooves allows knotless fixation of SERASYNTH® LOC. SERASYNTH® SLIM SLING is a reinforced suture consisting of several threads. The underlying suture material complies with the essential requirements of the Medical Device Directive 93/42/EEC and the European Pharmacopoeia (EP) monograph "Sutures, sterile synthetic absorbable (Fila resorbilia synthetica monofilamenta sterilia)", although SERASYNTH® LOC deviates from the pharmacopoeia in that it is the linear and not the knot tensile strength that forms the basis for suture stability.

INDICATIONS

SERASYNTH® is indicated for use in soft-tissue approximation and/or for ligation in cases where long-term (up to six weeks) absorbable suture material is indicated, including use in cardiovascular surgery, ophthalmic procedures and microsurgery.


SERASYNTH® LOC / SERASYNTH® SLIMSILING is intended for soft-tissue adaptation when long-term (up to six weeks) absorbable sutures are indicated.

USE

When selecting and using the suture the condition of the patient, the surgeon's experience, the surgical procedure, and the type and size of the wound should be taken into account. Knot security of the suture material requires the estab-

lished standard knotting techniques for synthetic monofilament sutures. To avoid damage to the needle while suturing it is recommended that the needle is grasped towards the needle base in an area between the end of the first third and the middle of the needle.

Endosuture SERASYNTH®: due to the material characteristics of the suture base material PDO, not all the conventional knots are suitable for the SERASYNTH® Endosuture. SERAG-WIESSNER has developed a novel knot which is marked by its ease of tying and good knot security. Please proceed as shown in the diagram:

ACTION

The purpose of the surgical suture is to create tissue-to-tissue approximation during the treatment of wounds. The ligature, or tie, as a special form of suture serves to occlude lumens.

The polymer PDO is degraded in the tissue by means of hydrolysis to 2-hydroxyethoxy acetic acid which is subsequently metabolised by the organism. This absorption expresses itself as a loss of tensile strength which is later accompanied by loss of mass. After 2 weeks tensile strength is reduced to approximately 70 – 80 % and after 4 weeks to 50 – 70 % of its original amount. Absorption is complete after about 30 weeks.

CONTRAINDICATIONS

SERASYNTH® / SERASYNTH® LOC / SERASYNTH® SLIMSILING should not be used where lasting stability of the suture material is required, as for instance when securing synthetic implants (heart valves, etc.).

WARNINGS / PRECAUTIONS / INTERACTIONS

The demands placed on the suture material and the required techniques will vary, depending on the area of application. The user should therefore be familiar with the various surgical techniques and procedures. Its in-vivo properties should be considered when selecting the suture material. The use should be carefully considered in patients with delayed wound healing. As with any foreign body there is a risk of calculus formation with prolonged contact of SERASYNTH® / SERASYNTH® LOC / SERASYNTH® SLIMSLING with salt solutions (e.g. in the biliary or urinary tracts). Sutures under tension or those needing additional support may require the supplemental use of non-absorbable sutures or suitable immobilisation procedures. In intracutaneous sutures the suture material – in particular the knots – should be placed as deep as possible. To ensure that the anchor fixation system functions properly, SERASYNTH® LOC should be handled carefully and the suture finished off with a sufficiently long backstitch, for example. There is still only limited experience of use in the nervous and cardiovascular systems, for intestinal anastomoses, and the closure of fascia. The surgical instruments should be handled as carefully as possible to avoid damage to the suture. Bending the needles or grasping them outside the recommended area may result in needle breakage. The handling and disposal of needles should be carried out with particular care due to the risk of contamination following needle stick injuries.

ADVERSE REACTIONS

Even with proper use the following adverse effects may occur:
Transient localised irritation, transient inflammatory foreign-body reactions, calculus formation with prolonged contact with salt solutions.

STERILITY

SERASYNTH® / SERASYNTH® LOC / SERASYNTH® SLIMSLING is supplied sterile (sterilization procedure: ethylene oxide gas). Do not resterilize,

critical changes in suitability for use possible (more details available on request)! Do not use damaged packages! Discard opened unused packs!

STORAGE CONDITIONS

No special storage conditions are required. Do not use after the expiry date.

DISPOSAL INFORMATION

Dispose of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

PRESENTATION

The products are available in various strengths and lengths attached to various needle types, and unneeded. The sutures and their combinations are supplied sterile in single and multiple packs (multipack). Please see catalogue for precise details.

EXPLANATION OF SYMBOLS ON THE PACKAGE

 1014 CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

Surgical suture, absorbable, monofilament, dyed

Polydioxanone

Antireflective needle

Detachable needle

Monodirectional anchor fixation system

Needle identification symbols and abbreviations are explained in the catalogue.

Do not use if package
is damaged

Do not re-use

Sterilized with ethylene oxide

MATERIAUX DE SUTURE en POLYDIOXANONE

- SERASYNTH®
- SERASYNTH® LOC
- SERASYNTH® SLIMSLING

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SERAG-WIESSNER GmbH & Co. KG