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● SERASYNTH® LOC

WARNINGS AND PRECAUTIONS

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified personnel, that is responsible to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: ethylene oxide).

STERILE EO

- The product is intended for single use only.

- Do not resterilize or re-use the product as critical changes in performance and safety are possible.

- Do not use after expiration date.

- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or unintentionally opened prior to use.

- Use the product only for the intended indication.

DESCRIPTION AND INTENDED PURPOSE

SERASYNTH® LOC is a synthetic, absorbable, monofilament surgical suture material. It is made of the polymer poly-p-dioxanone (PDO). It is available dyed with D&C violet No. 2, Colour Index No. 60725. SERASYNTH® LOC is available in various gauge sizes and lengths. It is supplied with various types of needles made of surgical steel. An anchor fixation system based on monodirectional barbs allows knotless fixation of SERASYNTH® LOC.

The underlying suture material complies with the current edition of the European Pharmacopoeia (EP) monograph for "Sutures, sterile synthetic absorbable (Fila resorbilia synthetica monofilamenta sterilia)", although SERASYNTH® LOC deviates from the pharmacopoeia in that it is the linear and not the knot tensile strength that forms the basis for suture stability.

SERASYNTH® LOC is intended for soft-tissue adaptation when long-term (up to six weeks) absorbable sutures are indicated.

INSTRUCTIONS FOR USE

In the choice and use of suture material, due account should be taken of the condition of the patient, the experience of the qualified personnel, the surgical procedure, and the type and size of the wound. In order to prevent needle damage during suturing it is

recommended that the needle be grasped at an area between the end of its first third (as seen from its end) and its midpoint.

MODE OF ACTION

The function of a medical suture in wound care is to attach tissue to tissue. The polymer PDO is degraded in the tissue by means of hydrolysis to 2-hydroxyethoxy acetic acid which is subsequently metabolised by the organism. This absorption expresses itself as a loss of tensile strength which is later accompanied by loss of mass. At least 50% of the initial tensile strength remains after 28 days (standardised in-vitro tests). Absorption is expected to be completed after about 182 days.

CONTRAINDICATIONS

SERASYNTH® LOC must not be used where long-term stability of the suture material is required. SERASYNTH® LOC is not indicated for use at the heart, in the circulatory or nervous system, in ophthalmology and microsurgery. The device must not be used in patients with allergies to any of the product components.

RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS

The requirements of the suture material and the techniques required vary with the indication. In the choice of suture material, account should be taken of the characteristics of the suture material concerned. The use of the device should be carefully considered in patients with delayed wound healing. Considerations should be taken regarding the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur. The use of absorbable sutures may be inappropriate in elderly, malnourished or debilitated patients. Like all foreign bodies, the device can give rise to calculus formation when in prolonged contact with salt-containing solutions (e.g. in the biliary or urinary tract). Additional use of nonabsorbable suture materials or suitable immobilization techniques may be required with sutures that are under tension or require greater support. In intracutaneous sutures the suture material should be placed as deep as possible. Skin sutures which remain in place longer than 7 days may cause local irritation and should be snipped off or removed. To ensure that the anchor fixation system functions properly, SERASYNTH® LOC should be handled carefully and the suture finished off with a sufficiently long backstitch, for example. Surgical instruments should be handled with care in order to avoid damage to the thread. Bending of needles or grasping them outside of the recommended area can result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled with particular care.

ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of the device: local irritation, inflammatory reactions to the foreign body, formation of foreign-body giant cells, suture fistulas or granulomas.

Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

STORAGE CONDITIONS

The storage conditions apply as indicated on the label.

DISPOSAL INFORMATION

Disposal of unused products or waste material with particular care due to the potential

risk of contamination and/or injury. In all other cases observe the specific national regulations.

MEANING OF THE SYMBOLS

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG