# SERATAN®

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## SUTURE MATERIAL MADE FROM TITANISED POLYAMIDE

### DESCRIPTION

SERATAN® is a synthetically manufactured, monofilament, non-absorbable, sterile surgical suture material. The basic suture material consists of polyamide and non-alloyed, oxidised titanium. Our polyamide suture materials are supplied as combinations with stainless steel needles. They are also supplied in different variants with specific accessories for intracutaneous sutures. They comply with the essential requirements of Directive 93/42/EEC concerning medical devices and the harmonised standard of the European Pharmacopoeia "Sutures, sterile, non-absorbable (Fila non resorbilia sterilia)".

### INDICATIONS

SERATAN® is indicated for soft-tissue adaptation and is, in principle, suitable for all wound sutures, in which non-absorbable thread material is indicated and no long-term stability of the suture material is required. Its pliability, low tendency to swell and the absence of wicking make SERATAN® an ideal suture material for all types of skin sutures. SERATAN® can also be used for marking purposes.

### USE

When selecting and using suture material, the patient's condition, the physician's experience, the surgical technique used and the wound size should be taken into consideration. Standard knot-tying techniques should be used to ensure knot-holding security. In order to prevent needle damage during suturing, it is recommended that the needle is grasped at a point between the end of its first third (as seen from its end) and the midpoint of the needle.

### ACTION

In wound care, the function of a medical suture is to attach tissue to tissue or tissue to allogeneic material.

Titanium is very well tolerated by the body and, as a material of choice in medicine, its use is favoured in many implants. The otherwise obligate tissue reaction to foreign bodies is much less severe if titanised

polyamide sutures are used. In the body, SERATAN® is very rapidly colonised by connective tissue cells and should therefore be removed at an early stage, normally approx 8 days, although this does also depend on its application site in the body (extremities, head, face etc.). It is non-absorbable, but can nevertheless lose its tensile strength over a long period and eventually fragment. As a result of its monofilament thread structure, SERATAN® can also be used in tissue at risk of infection.

### CONTRAINDICATIONS

SERATAN® is not indicated for use in the central vascular or central nervous systems. Due to the loss of tensile strength that occurs in the body, polyamide threads are contraindicated in cases where long-term stability of the suture material is required.

### WARNINGS / PRECAUTIONS / INTERACTIONS

The requirements on the suture material and the necessary techniques vary according to the indication. The user should therefore be experienced in the surgical techniques. When selecting the suture material, its in-vivo properties should be taken into account. For sutures placed in infected tissues, the use of appropriate concomitant measures is recommended. Surgical instruments should be handled carefully, in order to avoid damage to threads. This is generally of importance with monofilament suture material. Needles should be handled and disposed of with particular caution, due to the risk of contamination as a result of needle prick injuries.

### ADVERSE REACTIONS

Following correct use of SERATAN®, the following adverse reactions can occur: transient local irritation.

### STERILITY

The materials are delivered sterile (sterilization method: ethylene oxide). Do not resterilize, critical changes in suitability for use possible (more details available on request)!

Packs that have been opened for use but then not used are to be discarded!


### STORAGE CONDITIONS

The suture material should be stored in a clean dry place at room temperature. Do not use after the expiry date.

## PRESENTATION

SERATAN® is supplied in different gauge sizes and lengths in the form of needle-thread combinations with various needle types. The combinations are available in individual and multiple packs (Multi-pack) in sterile packaging and, in the case of the intracutaneous suture, with specific accessories. More precise details are given in the catalogue.

## EXPLANATION OF SYMBOLS ON THE PACKAGE

1014 CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

Surgical suture, nonabsorbable, monofilament, coated, undyed

Polyamide/titanium

Antireflective needle

Detachable needle

Needle identification symbols and abbreviations are explained in the catalogue.

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SERAG-WIESSNER GmbH & Co. KG