SERATAN®**de GEBRAUCHSANWEISUNG****NAHTMATERIAL AUS TITANISIERTEM POLYAMID****en INSTRUCTIONS FOR USE****SUTURE MATERIAL MADE FROM TITANISED POLYAMIDE****fr MODE D'EMPLOI****MATERIAU DE SUTURE EN POLYAMIDE TITANISE****it ISTRUZIONI PER L'USO****MATERIALE DA SUTURA IN POLIAMMIDE TITANIZZATA****es INSTRUCCIONES DE USO****MATERIAL DE SUTURA DE POLIAMIDA TITANIZADA****nl GEBRUIKSAANWIJZING****HECHTMATERIAAL UIT GETITANISEERD POLYAMIDE****pl INSTRUKCJA OBSŁUGI****NICI CHIRURGICZNE Z TYTANOWANEGO POLIAMIDU****cs NÁVOD K POUŽITÍ****ŠICÍ MATERIÁL Z TITANIZOVANÉHO POLYAMIDU****tr KULLANIM TALİMATI****TITANIZE EDILMIŞ POLİAMİDDEN DİKİŞ MATERYALİ****lt NAUDOJIMO INSTRUKCIJA****SIŪLAI IŠ TITANU DENGTO POLIAMIDO****lv LIETOŠANAS INSTRUKCIJA****ĶĪRURĢISKIE ŠŪJAMIE MATERIĀLI NO TITANISIERTEM POLYAMID**

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SUTURE MATERIAL MADE FROM TITANISED POLYAMIDE

WARNINGS AND PRECAUTIONS

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified and trained personnel. It is the physician's responsibility to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: ethylene oxide). 
- The product is intended for single use only. 
- Do not re-sterilize or reuse the product as critical changes in performance and safety are possible. 
- Do not use after expiration date. 
- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or accidentally opened prior to use. 
- Use the product only for the intended indication.
- Please contact the manufacturer for further information.

DESCRIPTION AND INTENDED PURPOSE

SERATAN[®] is a monofilament, non-absorbable, surgical suture material. The basic suture material consists of polyamide and non-alloyed, oxidised titanium. SERATAN[®] is supplied dyed blue with Phthalocyanine blue, Colour Index No. 74160 or undyed. Our suture materials are available in various gauge sizes and lengths. They are supplied in combinations with various types of needles made of surgical steel. They comply with the current edition of the European Pharmacopoeia (EP) monograph "Sutures, sterile non-absorbable" (fila non resorbilia sterilia)".

SERATAN[®] is indicated for soft-tissue adaptation and is, in principle, suitable for all wound sutures, in which non-absorbable thread material is indicated and no long-term stability of the suture material is required. Threads left in tissues for a long time can lose tensile strength and eventually can fragment. SERATAN[®] is mainly used for skin sutures. SERATAN[®] can also be used for marking purposes.

Titanium is well tolerated and, as a material of choice in medicine, its use is favoured in many implants. The tissue reaction to foreign bodies is less severe if titanised polyamide sutures are used.

INSTRUCTIONS FOR USE

In the choice and use of suture material, due account should be taken of the condition of the patient, the experience of the physician, the surgical technique, and the type and size of the wound. Standard knot-tying techniques should be used to ensure knot-holding security.

In order to prevent needle damage during suturing it is recommended that the needle be grasped at an area between the end of its first third (as seen from its end) and its

midpoint.

ACTION

The function of a medical suture in wound care is to attach tissue to tissue or tissue to foreign material. SERATAN® is very rapidly colonised by connective tissue cells and should therefore be removed at an early stage, usually within approx. 8 days. This does also depend on its application site in the body (extremities, head, face etc.).

CONTRAINDICATIONS

SERATAN® is not indicated for use in the central vascular or central nervous systems. Due to the loss of tensile strength that occurs in the body, polyamide threads are contraindicated in cases where long-term stability of the suture material is required. The device must not be used in patients with allergies to any of the product components.

RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS

The requirements of the suture material and the techniques required vary with the indication. In the choice of suture material, account should be taken of the characteristics of the suture material concerned. In view of the high elasticity of thicker titanised monofilament polyamide threads, it is recommended that one to two additional knots be applied in order to achieve knot-holding security.

As a result of its monofilament thread structure, SERATAN® can also be used in tissue at risk of infection. Use of suitable ancillary measures is recommended with sutures in infected tissue. Surgical instruments should be handled with care in order to avoid damage to the thread. This is generally of importance with monofilament suture material. Bending of needles or grasping them outside of the recommended area can result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled with particular care.

ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of the device: local irritation, inflammatory reactions to the foreign body; formation of foreign-body giant cells, suture fistulas or granulomas, allergic reactions to a component of the product.

Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

STORAGE CONDITIONS

The storage conditions apply as indicated on the label.

DISPOSAL INFORMATION

Disposal of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

EXPLANATION OF SYMBOLS

Symbols are explained at the end of the instructions for use.

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**SERAG-WIESSNER GmbH & Co. KG**